

For Research Use Only

RPUSD4 Polyclonal antibody

Catalog Number: 27843-1-AP



Basic Information

Catalog Number: 27843-1-AP	GenBank Accession Number: BC014131	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 84881	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q96CM3	
Isotype: IgG	Full Name: RNA pseudouridylylate synthase domain containing 4	
Immunogen Catalog Number: AG27225	Observed MW: 39-42 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, rat liver tissue
Species Specificity: Human, Mouse, Rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27843-1-AP (RPUSD4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.