

For Research Use Only

Maturin Polyclonal antibody

Catalog Number: 27936-1-AP



Basic Information

Catalog Number:

27936-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19932

GenBank Accession Number:

BC156839

GeneID (NCBI):

222166

UNIPROT ID:

Q8N3F0

Full Name:

chromosome 7 open reading frame 41

Calculated MW:

131 aa, 15 kDa

Observed MW:

14-17 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, pig, rat

Positive Controls:

WB : COLO 320 cells, pig brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27936-1-AP (Maturin antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 27936-1-AP (Maturin antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.