

For Research Use Only

Maturin Polyclonal antibody

Catalog Number: 27936-1-AP



Basic Information

Catalog Number: 27936-1-AP	GenBank Accession Number: BC156839	Purification Method: Antigen affinity purification
Size: 150 µg/ml	GeneID (NCBI): 222166	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8N3F0	
Isotype: IgG	Full Name: chromosome 7 open reading frame 41	
Immunogen Catalog Number: AG19932	Calculated MW: 131 aa, 15 kDa	
	Observed MW: 14-17 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COLO 320 cells, pig brain tissue
Species Specificity: Human, mouse, pig, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27936-1-AP (Maturin antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 27936-1-AP (Maturin antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.