

For Research Use Only

# ZBTB24 Polyclonal antibody

Catalog Number: 27975-1-AP

1 Publications



## Basic Information

Catalog Number:

27975-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24458

GenBank Accession Number:

BC036731

GeneID (NCBI):

9841

UNIPROT ID:

O43167

Full Name:

zinc finger and BTB domain containing 24

Calculated MW:

697 aa, 78 kDa

Observed MW:

80-85 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : HeLa cells, Jurkat cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Noriko Miyake	31839203	Am J Hum Genet	WB,IP

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27975-1-AP (ZBTB24 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.