

For Research Use Only

RHNO1 Polyclonal antibody, PBS Only

Catalog Number: 27977-1-PBS



Basic Information

Catalog Number:

27977-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14789

GenBank Accession Number:

BC005106

GeneID (NCBI):

83695

Full Name:

chromosome 12 open reading frame

32

Calculated MW:

238 aa, 27 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

RHNO1 (Rad9-Hus1-Rad1 Interacting Nuclear Orphan 1) is a protein associated with the DNA replication stress (DRS) response and DNA repair. It interacts with the 9-1-1 checkpoint clamp and TopBP1 to activate the ATR/Chk1 signaling pathway, thereby maintaining genomic integrity. Additionally, RHNO1 has been identified as a key facilitator of theta-mediated end joining (TMEJ), a DNA repair mechanism related to cancer progression and chemotherapy resistance. In cancer research, RHNO1 is considered a potential therapeutic target due to its high expression in tumor cells and its regulatory role in cell proliferation and survival.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-RHNO1 (IP:27977-1-AP, 4ug; Detection:27977-1-AP 1:500) with mouse thymus tissue lysate 1320 ug. This data was developed using the same antibody clone with 27977-1-PBS in a different storage buffer formulation.