

For Research Use Only

ABLIM3 Polyclonal antibody

Catalog Number: 27981-1-AP

1 Publications



Basic Information

Catalog Number:

27981-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27578

GenBank Accession Number:

BC001665

GeneID (NCBI):

22885

UNIPROT ID:

O94929

Full Name:

actin binding LIM protein family,
member 3

Calculated MW:

78 kDa

Observed MW:

78 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse

Positive Controls:

WB : NIH/3T3 cells, mouse brain tissue, HepG2 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Joleen S Cheah	33909446	Mol Biol Cell	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27981-1-AP (ABLIM3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.