

For Research Use Only

Tissue Factor/CD142 Polyclonal antibody, PBS Only

Catalog Number: 28005-1-PBS



Basic Information

Catalog Number:

28005-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27553

GenBank Accession Number:

BC011029

GeneID (NCBI):

2152

UNIPROT ID:

P13726

Full Name:

coagulation factor III (thromboplastin, tissue factor)

Calculated MW:

295 aa, 33 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Normal blood coagulation is a complex process, involving a cascade of activation of different plasma proteins, ultimately resulting in the formation of a clot, called fibrin. Tissue Factor (TF), also named Coagulation factor III or F3, is the primary initiator of the blood coagulation cascade and plays an essential role in hemostasis (PMID: 26877187; PMID: 33796108).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28005-1-AP (Tissue Factor antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 28005-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 28005-1-AP (Tissue Factor antibody) at dilution of 1:100 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 28005-1-PBS in a different storage buffer formulation.