

For Research Use Only

NDUFAF7 Polyclonal antibody

Catalog Number: 28025-1-AP

1 Publications



Basic Information

Catalog Number:

28025-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27590

GenBank Accession Number:

BC004548

GeneID (NCBI):

55471

UNIPROT ID:

Q7L592

Full Name:

chromosome 2 open reading frame 56

Calculated MW:

441 aa, 49 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Mouse, Rat, human

Cited Species:

mouse

Positive Controls:

WB : mouse heart tissue, rat heart tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ahmed M Elmansi	37273159	Geroscience	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28025-1-AP (C2orf56 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.