

For Research Use Only

SMG8 Polyclonal antibody

Catalog Number: 28040-1-AP



Basic Information

Catalog Number:

28040-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27820

GenBank Accession Number:

BC020957

GeneID (NCBI):

55181

UNIPROT ID:

Q8ND04

Full Name:

chromosome 17 open reading frame

71

Observed MW:

110-120 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28040-1-AP (SMG8 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.