

For Research Use Only

PPARD Polyclonal antibody

Catalog Number: 28053-1-AP

4 Publications



Basic Information

Catalog Number:

28053-1-AP

Concentration:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27817

GenBank Accession Number:

BC002715

GeneID (NCBI):

5467

UNIPROT ID:

Q03181

Full Name:

peroxisome proliferator-activated
receptor delta

Calculated MW:

50 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

Human

Cited Species:

human, bovine

Positive Controls:

WB : human placenta tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jing Qian	39993702	Eur J Pharmacol	WB
Yilei Guo	39632616	J Cell Biochem	WB, IHC
Shanshan Yang	39765775	Antioxidants (Basel)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 28053-1-AP (PPARD antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.