

For Research Use Only

CCDC45 Polyclonal antibody

Catalog Number: 28062-1-AP



Basic Information

Catalog Number:

28062-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27821

GenBank Accession Number:

BC009518

GeneID (NCBI):

90799

UNIPROT ID:

Q96GE4

Full Name:

coiled-coil domain containing 45

Calculated MW:

821 aa, 95 kDa

Observed MW:

95-105 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Jurkat cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28062-1-AP (CCDC45 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.