

For Research Use Only

# PEAR1 Polyclonal antibody

Catalog Number: 28086-1-AP



## Basic Information

Catalog Number:

28086-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27978

GenBank Accession Number:

NM\_001080471

GeneID (NCBI):

375033

UNIPROT ID:

Q5VY43

Full Name:

platelet endothelial aggregation  
receptor 1

Calculated MW:

111 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Human peripheral blood platelets tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 28086-1-AP (PEAR1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.