

For Research Use Only

# ASIC3 Polyclonal antibody

Catalog Number: 28093-1-AP



## Basic Information

Catalog Number:

28093-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26092

GenBank Accession Number:

NM\_004769

GeneID (NCBI):

9311

UNIPROT ID:

Q9UHC3

Full Name:

amiloride-sensitive cation channel 3

Calculated MW:

59 kDa

Observed MW:

59-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Mouse kidney tissue lysate was subjected to SDS PAGE followed by western blot with 28093-1-AP (ASIC3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.