

For Research Use Only

ZC3H14 Polyclonal antibody, PBS Only

Catalog Number: 28190-1-PBS



Basic Information

Catalog Number:

28190-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28182

GenBank Accession Number:

NM_001160103

GeneID (NCBI):

79882

UNIPROT ID:

Q6PJT7

Full Name:

zinc finger CCCH-type containing 14

Calculated MW:

83 aa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 28190-1-AP (ZC3H14 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 28190-1-PBS in a different storage buffer formulation.



IP result of anti-ZC3H14 (IP:28190-1-AP, 4ug; Detection:28190-1-AP 1:4000) with HEK-293T cells lysate 2200 ug. This data was developed using the same antibody clone with 28190-1-PBS in a different storage buffer formulation.