

SFRP1 Polyclonal antibody

Catalog Number: 28271-1-AP

Basic Information

Catalog Number:

28271-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28072

GenBank Accession Number:

BC036503

GeneID (NCBI):

6422

UNIPROT ID:

Q8N474

Full Name:

secreted frizzled-related protein 1

Calculated MW:

314 aa, 35 kDa

Observed MW:

30-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : human testis tissue,

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 28271-1-AP (SFRP1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.