

For Research Use Only

Phospho-DLK (Ser272) Polyclonal antibody



Catalog Number: 28279-1-AP

Basic Information

Catalog Number: 28279-1-AP	GenBank Accession Number: NM_001193511	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 7786	
Source: Rabbit	UNIPROT ID: Q12852	
Isotype: IgG	Full Name: mitogen-activated protein kinase kinase kinase 12	
	Calculated MW: 93 kda	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data