

For Research Use Only

# PDS5B Polyclonal antibody

Catalog Number: 28318-1-AP



## Basic Information

Catalog Number:

28318-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27977

GenBank Accession Number:

NM\_015032

GeneID (NCBI):

23047

UNIPROT ID:

Q9NTI5

Full Name:

PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae)

Calculated MW:

165 aa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : PC-3 cells, SKOV-3 cells, mouse testis tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28318-1-AP (PDS5B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.