

For Research Use Only

# PDS5B Polyclonal antibody

Catalog Number: 28318-1-AP



## Basic Information

**Catalog Number:**

28318-1-AP

**Size:**

1000 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG27977

**GenBank Accession Number:**

NM\_015032

**GeneID (NCBI):**

23047

**UNIPROT ID:**

Q9NTI5

**Full Name:**

PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae)

**Calculated MW:**

165 aa

**Observed MW:**

150 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

WB 1:1000-1:4000

## Applications

**Tested Applications:**

WB, ELISA

**Species Specificity:**

Human, mouse

**Positive Controls:**

WB: PC-3 cells, SKOV-3 cells, mouse testis tissue

## Background Information

### Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28318-1-AP (PDS5B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.