

For Research Use Only

# ZC3H13 Polyclonal antibody

Catalog Number: 28430-1-AP



## Basic Information

Catalog Number:

28430-1-AP

Size:

330 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28909

GenBank Accession Number:

NM\_001076788

GeneID (NCBI):

23091

UNIPROT ID:

Q5T200

Full Name:

zinc finger CCCH-type containing 13

Calculated MW:

197 kDa

Observed MW:

196-240 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293T cells, HeLa cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28430-1-AP (ZC3H13 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.