

For Research Use Only

CCDC171 Polyclonal antibody

Catalog Number: 28433-1-AP



Basic Information

Catalog Number:

28433-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28926

GenBank Accession Number:

NM_173550

GeneID (NCBI):

203238

UNIPROT ID:

Q6TFL3

Full Name:

chromosome 9 open reading frame 93

Calculated MW:

153 kDa

Observed MW:

142 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat, human

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28433-1-AP (CCDC171 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.