

For Research Use Only

CCDC171 Polyclonal antibody

Catalog Number: 28433-1-AP



Basic Information

Catalog Number: 28433-1-AP	GenBank Accession Number: NM_173550	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 203238	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q6TFL3	
Isotype: IgG	Full Name: chromosome 9 open reading frame 93	
Immunogen Catalog Number: AG28926	Calculated MW: 153 kDa	
	Observed MW: 142 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Species Specificity: mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28433-1-AP (CCDC171 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.