

For Research Use Only

NOC2L Polyclonal antibody

Catalog Number: 28509-1-AP



Basic Information

Catalog Number: 28509-1-AP	GenBank Accession Number: BC003555	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 26155	Recommended Dilutions: WB 1:2000-1:12000 IHC 1:500-1:2000 IF 1:50-1:500
Source: Rabbit	UNIPROT ID: Q9Y3T9	
Isotype: IgG	Full Name: nucleolar complex associated 2 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG26952	Calculated MW: 749 aa, 85 kDa	
	Observed MW: 85 kDa	

Applications

Tested Applications:
IF/ICC, IHC, WB, ELISA

Species Specificity:
Human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: U2OS cells, HEK-293T cells, HuH-7 cells, mouse cerebellum tissue, mouse stomach tissue

IHC: human prostate cancer tissue,

IF: U2OS cells,

Background Information

NOC2L, also known as NIR, is a transcriptional corepressor with inhibitor of histone acetyltransferase (INHAT) activity (PMID: 16322561). It prevents acetylation of all core histones by the histone acetyltransferase EP300/p300. NOC2L interacts with p53, inhibits p53-activated gene expression, and regulates p53-dependent apoptosis (PMID: 16322561). NIR binds to TAp63 and acts as a repressor of TAp63-mediated transactivation of the p21Cip1 gene (PMID: 20123734). NOC2L is ubiquitously expressed throughout embryonic development and adulthood (PMID: 16322561). NOC2L has been reported as a key regulator of skin development. Epidermis-specific ablation of NOC2L impairs epidermal stratification and barrier function (PMID: 29490983).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 28509-1-AP (NOC2L antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 28509-1-AP (NOC2L antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 28509-1-AP (NOC2L antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using 28509-1-AP (NOC2L antibody), at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Phalloidine stains the F-actin (red).