

For Research Use Only

PLAG1 Polyclonal antibody

Catalog Number: 28651-1-AP



Basic Information

Catalog Number:

28651-1-AP

Concentration:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29010

GenBank Accession Number:

BC075047

GeneID (NCBI):

5324

UNIPROT ID:

Q6DJT9

Full Name:

pleiomorphic adenoma gene 1

Calculated MW:

500 aa, 56 kDa

Observed MW:

56-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28651-1-AP (PLAG1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.