

For Research Use Only

CCDC60 Polyclonal antibody

Catalog Number: 28652-1-AP



Basic Information

Catalog Number:

28652-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28109

GenBank Accession Number:

BC040553

GeneID (NCBI):

160777

UNIPROT ID:

Q8IWA6

Full Name:

coiled-coil domain containing 60

Calculated MW:

550 aa, 63 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HepG2 cells,

Background Information

Coiled-coil domain containing 60 (CCDC60) is a member of the CCDC family, which participates in the progression of many types of cancer.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 28652-1-AP (CCDC60 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.