

For Research Use Only

TRABD2B Polyclonal antibody

Catalog Number: 28682-1-AP



Basic Information

Catalog Number:

28682-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30248

GenBank Accession Number:

NM_001194986

GeneID (NCBI):

388630

UNIPROT ID:

A6NFA1

Full Name:

UPF0632 protein A

Calculated MW:

57 kDa

Observed MW:

65-79 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : DU 145 cells, HT-29 cells, HeLa cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28682-1-AP (TRABD2B antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.