

For Research Use Only

# PCAF Polyclonal antibody

Catalog Number: 28770-1-AP **5 Publications**



## Basic Information

Catalog Number:

28770-1-AP

Concentration:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30238

GenBank Accession Number:

BC060823

GeneID (NCBI):

8850

UNIPROT ID:

Q92831

Full Name:

K(lysine) acetyltransferase 2B

Calculated MW:

93 kDa

Observed MW:

95-110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

Human

Cited Species:

human, mouse

Positive Controls:

WB : A431 cells, HeLa cells

## Background Information

### Notable Publications

| Author       | Pubmed ID | Journal        | Application |
|--------------|-----------|----------------|-------------|
| Lipeng Pei   | 39901144  | Mol Med        | WB          |
| Enyuan Zhang | 39012622  | J Proteome Res | WB          |
| Yihong Li    | 38821061  | Cancer Cell    | WB          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28770-1-AP (KAT2B antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.