

For Research Use Only

AMD1 Polyclonal antibody

Catalog Number: 29108-1-AP



Basic Information

Catalog Number:

29108-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30751

GenBank Accession Number:

BC000171

GeneID (NCBI):

262

UNIPROT ID:

P17707

Full Name:

adenosylmethionine decarboxylase 1

Calculated MW:

38 kDa, 32 kDa

Observed MW:

30-40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : A549 cells, DU 145 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29108-1-AP (AMD1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.