

For Research Use Only

# DHX58 Polyclonal antibody

Catalog Number: 29114-1-AP



## Basic Information

Catalog Number:

29114-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30760

GenBank Accession Number:

BC014949

GeneID (NCBI):

79132

UNIPROT ID:

Q96C10

Full Name:

DEXH (Asp-Glu-X-His) box polypeptide 58

Calculated MW:

678 aa, 77 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

