

For Research Use Only

ALG9 Polyclonal antibody, PBS Only

Catalog Number: 29218-1-PBS



Basic Information

Catalog Number:

29218-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30489

GenBank Accession Number:

BC009255

GeneID (NCBI):

79796

UNIPROT ID:

Q9H6U8

Full Name:

asparagine-linked glycosylation 9,
alpha-1,2-mannosyltransferase
homolog (S. cerevisiae)

Calculated MW:

70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

ALG9 encodes an endoplasmic reticulum enzyme that builds N-glycans, the third ER protein-encoding polycystic disease gene after GANAB and DNAJB11. Autosomal recessive loss of ALG9 results in a severe congenital disorder of glycosylation (CDG) with a multiorgan phenotype that includes kidney cysts (PMID: 31395617). Western blot analysis detected ALG9 at an apparent molecular mass of 70 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29218-1-AP (ALG9 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 29218-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 29218-1-AP (ALG9 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 29218-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 29218-1-AP (ALG9 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 29218-1-PBS in a different storage buffer formulation.