

For Research Use Only

# SARS-COV-2 NSP7 Polyclonal antibody

Catalog Number: 29263-1-AP



## Basic Information

Catalog Number:

29263-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30933

GenBank Accession Number:

YP\_009742614

GeneID (NCBI):

43740578

Full Name:

ORF1ab

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

Virus

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



SARS-CoV-2 Envelope protein Antibody (29263-1-AP) tested by ELISA. SARS-COV-2 NSP7 protein was coated onto microtiter plates at 1 µg/well and then incubated with a dilution series of SARS-COV-2 NSP7 protein Antibody (29263-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.