

For Research Use Only

ZNF664 Polyclonal antibody

Catalog Number: 29317-1-AP



Basic Information

Catalog Number:

29317-1-AP

Size:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30992

GenBank Accession Number:

BC101530

GeneID (NCBI):

144348

UNIPROT ID:

Q8N3J9

Full Name:

zinc finger protein 664

Calculated MW:

261 aa, 30 kDa

Observed MW:

30-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

IP : mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ZNF664 (IP:29317-1-AP, 4ug;
Detection:29317-1-AP 1:500) with mouse liver
tissue lysate 2400 ug.