

For Research Use Only

C19orf72 Polyclonal antibody

Catalog Number: 29362-1-AP



Basic Information

Catalog Number:

29362-1-AP

Size:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30318

GenBank Accession Number:

BC013280

GeneID (NCBI):

90379

UNIPROT ID:

Q66K64

Full Name:

chromosome 19 open reading frame

72

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : DU 145 cells, K-562 cells, PC-3 cells

Background Information

C19orf72, also called DCAF15 (DDB1- and CUL4-associated factor 15), is the substrate-recognition component of the DCX complex, which is a cullin-4-RING E3 ubiquitin-protein ligase complex that mediates ubiquitination and degradation of target proteins (PMID:16949367; 31452512).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29362-1-AP (C19orf72 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.