

For Research Use Only

# SOAT1 Polyclonal antibody

Catalog Number: 29369-1-AP



## Basic Information

Catalog Number:

29369-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30945

GenBank Accession Number:

BC028940

GeneID (NCBI):

6646

UNIPROT ID:

P35610

Full Name:

sterol O-acyltransferase 1

Calculated MW:

550 aa, 65 kDa

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, SW480 cells, mouse liver tissue, rat liver tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29369-1-AP (SOAT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.