

For Research Use Only

P54 of ASFV Polyclonal antibody

Catalog Number: 29511-1-AP



Basic Information

Catalog Number:

29511-1-AP

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31444

GenBank Accession Number:

AAA65354

GeneID (NCBI):

22220355

UNIPROT ID:

Q65194

Full Name:

Envelope protein p54

Calculated MW:

20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

ELISA 1:556544-1:2226176

Applications

Tested Applications:

ELISA

Species Specificity:

human

Positive Controls:

ELISA : FusionProtein,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



P54 of ASFV Antibody (29511-1-AP) tested by ELISA. P54 of ASFV(1-65aa) were coated onto microtiter plates at 0.15 µg/well and then incubated with a dilution series of P54 of ASFV Antibody (29511-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.