

For Research Use Only

P30 of ASFV Polyclonal antibody

Catalog Number: 29513-1-AP



Basic Information

Catalog Number: 29513-1-AP	GenBank Accession Number: M96354	Purification Method: Antigen affinity purification
Size: 1200 µg/ml	GeneID (NCBI): 22220322	Recommended Dilutions: ELISA 1:256000-1:1024000
Source: Rabbit	UNIPROT ID: P34204	
Isotype: IgG	Full Name: Phosphoprotein p30	
Immunogen Catalog Number: AG31446	Calculated MW: 30 kDa	

Applications

Tested Applications: ELISA	Positive Controls: ELISA : FusionProtein,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



P30 of ASFV Antibody (29513-1-AP) tested by ELISA. P30 of ASFV(1-204aa) were coated onto microtiter plates at 0.15 µg/well and then incubated with a dilution series of P30 of ASFV Antibody (29513-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.