

For Research Use Only

# P30 of ASFV Polyclonal antibody

Catalog Number: 29513-1-AP



## Basic Information

Catalog Number:

29513-1-AP

Size:

1200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31446

GenBank Accession Number:

M96354

GeneID (NCBI):

22220322

UNIPROT ID:

P34204

Full Name:

Phosphoprotein p30

Calculated MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

ELISA 1:256000-1:1024000

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

Positive Controls:

ELISA : FusionProtein,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



P30 of ASFV Antibody (29513-1-AP) tested by ELISA. P30 of ASFV(1-204aa ) were coated onto microtiter plates at 0.15 µg/well and then incubated with a dilution series of P30 of ASFV Antibody (29513-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.