

For Research Use Only

SPINK5 Polyclonal antibody

Catalog Number: 29808-1-AP

1 Publications



Basic Information

Catalog Number:

29808-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30432

GenBank Accession Number:

NM_001127698

GeneID (NCBI):

11005

UNIPROT ID:

Q9NQ38

Full Name:

serine peptidase inhibitor, Kazal type 5

Calculated MW:

121 kDa

Observed MW:

120-130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB: COLO 320 cells,

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|---------------|-------------|
| Yaning Guan | 41080545 | Front Immunol | |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 29808-1-AP (SPINK5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.