

For Research Use Only

ESCO1 Polyclonal antibody

Catalog Number: 29821-1-AP



Basic Information

Catalog Number:

29821-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31189

GenBank Accession Number:

NM_052911

GeneID (NCBI):

114799

UNIPROT ID:

Q5FWF5

Full Name:

establishment of cohesion 1 homolog
1 (S. cerevisiae)

Calculated MW:

95 kDa

Observed MW:

95 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : A375 cells, COLO 320 cells, SGC-7901 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29821-1-AP (ESCO1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.