

For Research Use Only

ESCO1 Polyclonal antibody, PBS Only

Catalog Number:29821-1-PBS



Basic Information

Catalog Number:

29821-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31189

GenBank Accession Number:

NM_052911

GeneID (NCBI):

114799

UNIPROT ID:

Q5FWF5

Full Name:

establishment of cohesion 1 homolog
1 (S. cerevisiae)

Calculated MW:

95 kDa

Observed MW:

95 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29821-1-AP (ESCO 1 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 29821-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse small intestine tissue slide using 29821-1-AP (ESCO 1 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 29821-1-PBS in a different storage buffer formulation.