

For Research Use Only

TSR1 Polyclonal antibody

Catalog Number: 29957-1-AP



Basic Information

Catalog Number:

29957-1-AP

Size:

650 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31575

GenBank Accession Number:

NM_018128

GeneID (NCBI):

55720

UNIPROT ID:

Q2NL82

Full Name:

TSR1, 20S rRNA accumulation,
homolog (S. cerevisiae)

Calculated MW:

92KD

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

IF/ICC 1:200-1:800

Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells, L02 cells, U-251 cells

IF/ICC : U-251 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U-251 cells using TSR1 antibody (29957-1-AP) at dilution of 1:400 and Coralite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).

Various lysates were subjected to SDS PAGE followed by western blot with 29957-1-AP (TSR1 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.