

For Research Use Only

TSR1 Polyclonal antibody

Catalog Number: 29957-1-AP



Basic Information

Catalog Number:

29957-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31575

GenBank Accession Number:

NM_018128

GeneID (NCBI):

55720

UNIPROT ID:

Q2NL82

Full Name:

TSR1, 20S rRNA accumulation,
homolog (S. cerevisiae)

Calculated MW:

92KD

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293 cells, HeLa cells, L02 cells, U-251 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 29957-1-AP (TSR1 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.