

For Research Use Only

KIAA1797 Polyclonal antibody, PBS Only

Catalog Number: 30140-1-PBS

Featured Product



Basic Information

Catalog Number:

30140-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31233

GenBank Accession Number:

BC001246

GeneID (NCBI):

54914

UNIPROT ID:

Q5VW36

Full Name:

KIAA1797

Observed MW:

160-170 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

KIAA1797 is a novel component of the focal adhesion complex expressed in glial and neuronal cells with a tumour suppressor function in gliomas, which is also named 'focadhesin' (PMID: 22427331).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30140-1-AP (KIAA1797 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30140-1-PBS in a different storage buffer formulation.



WB result of KIAA1797 antibody (30140-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-KIAA1797 transfected HeLa cells. This data was developed using the same antibody clone with 30140-1-PBS in a different storage buffer formulation.