

For Research Use Only

IL-1 alpha Polyclonal antibody

Catalog Number:30216-1-AP



Basic Information

Catalog Number:

30216-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0263

GenBank Accession Number:

NM_010554

GeneID (NCBI):

16175

UNIPROT ID:

P01582

Full Name:

interleukin 1 alpha

Calculated MW:

31kd

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : LPS treated RAW 264.7 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



LPS treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 30216-1-AP (IL1-alpha antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.