

For Research Use Only

AKAP9 Polyclonal antibody, PBS Only

Catalog Number:30290-1-PBS



Basic Information

Catalog Number:

30290-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32425

GenBank Accession Number:

NM_005751.4

GeneID (NCBI):

10142

UNIPROT ID:

Q99996

Full Name:

A kinase (PRKA) anchor protein
(yotiao) 9

Calculated MW:

453 kDa

Observed MW:

450-500 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 30290-1-AP (AKAP9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30290-1-PBS in a different storage buffer formulation.



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 30290-1-AP (AKAP9 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30290-1-PBS in a different storage buffer formulation.