

For Research Use Only

# RB1 Polyclonal antibody

Catalog Number: 30374-1-AP



## Basic Information

Catalog Number:

30374-1-AP

Size:

1300 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC039060

GeneID (NCBI):

5925

UNIPROT ID:

P06400

Full Name:

retinoblastoma 1

Calculated MW:

928 aa, 106 kDa

Observed MW:

106-110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : Jurkat cells,

## Background Information

The RB1 gene, a member of the pocket family with p107 and p130, was the first known tumor suppressor. RB1 protein is targeted by the transforming proteins of the DNA tumor viruses like adenoviral E1A and is functionally inactivated in most human tumor cells owing to mutations of either the RB1 gene itself or its upstream regulators (PMID: 22787429).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30374-1-AP (RB1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.