

For Research Use Only

PCYOX1L Polyclonal antibody

Catalog Number: 30485-1-AP



Basic Information

Catalog Number:

30485-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33077

GenBank Accession Number:

BC000014

GeneID (NCBI):

78991

UNIPROT ID:

Q8NBM8

Full Name:

prenylcysteine oxidase 1 like

Calculated MW:

494 aa, 55 kDa

Observed MW:

55-57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : LNCaP cells, mouse cerebellum tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30485-1-AP (PCYOX1L antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 30485-1-AP (PCYOX1L antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.