

For Research Use Only

# SLX4IP Polyclonal antibody, PBS Only

Catalog Number:30497-1-PBS



## Basic Information

Catalog Number:

30497-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33034

GenBank Accession Number:

GeneID (NCBI):

128710

UNIPROT ID:

Q5VYV7

Full Name:

chromosome 20 open reading frame

94

Calculated MW:

46 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

SLX4IP is a nuclear adaptor protein that associates with the SLX4 DNA repair scaffold, contributing to DNA damage repair, replication stress responses, and telomere maintenance. By regulating nuclease-containing repair complexes and potentially supporting ALT-mediated telomere elongation, SLX4IP plays an important role in maintaining genomic stability. Dysregulation of SLX4IP has been implicated in cancer biology, particularly in tumors utilizing alternative telomere maintenance pathways.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30497-1-AP (SLX4IP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30497-1-PBS in a different storage buffer formulation.