

For Research Use Only

C2CD3 Polyclonal antibody

Catalog Number: 30541-1-AP



Basic Information

Catalog Number:

30541-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30423

GenBank Accession Number:

NM_001286577

GeneID (NCBI):

26005

UNIPROT ID:

Q4AC94

Full Name:

C2 calcium-dependent domain
containing 3

Calculated MW:

260 kDa

Observed MW:

260 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : mouse kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30541-1-AP (C2CD3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.