

For Research Use Only

CRKL Polyclonal antibody

Catalog Number: 30546-1-AP



Basic Information

Catalog Number:

30546-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30874

GenBank Accession Number:

NM_005207

GeneID (NCBI):

1399

UNIPROT ID:

P46109

Full Name:

v-crk sarcoma virus CT10 oncogene
homolog (avian)-like

Calculated MW:

34 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293T cells, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30546-1-AP (CRKL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.