

For Research Use Only

# LAPTM4A Polyclonal antibody, PBS Only

Catalog Number:30627-1-PBS



## Basic Information

Catalog Number:

30627-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6501

GenBank Accession Number:

BC000421

GeneID (NCBI):

9741

UNIPROT ID:

Q15012

Full Name:

lysosomal protein transmembrane 4

alpha

Calculated MW:

27 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IF-P, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

LAPTM4A is an intracellular membrane protein that regulates the function of organic cation transporter-2 (OCT2) through direct interaction. LAPTM4A can also form a complex with A4GALT in the Golgi and shows at about 70 kDa.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 30627-1-AP (LAPTM4A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 30627-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed mouse kidney tissue using LAPTM4A antibody (30627-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 30627-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed mouse kidney tissue using LAPTM4A antibody (30627-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 30627-1-PBS in a different storage buffer formulation.